

The RDS 6110 and RDS 3110 use the proprietary AeroClave process to decontaminate spaces.

The RDS 6110 can effectively treat spaces, including irregular shaped rooms, up to 1500 square feet such as gyms, locker rooms, classrooms and other institutional facilities.

The RDS 3110 is a ruggedized, self-contained, 48 lbs. portable decontamination system that can treat spaces up to 5,000 cubic feet like ambulances, cars/trucks, and small offices.

The RDS 6110 and RDS 3110 can be operated in either an aerosolized application mode for hands-free room decontamination or hand-applied mode. The disinfectant solution (Vital Oxide) is ready-to-use and produces no harmful by-products, does not contain volatile organic compounds, and is 100% biodegradable.

The AeroClave process, using the Vital Oxide disinfectant, is effective against Novel Coronavirus (2019-nCoV) based on CDC recommendations.

